



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)  
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,  
 68000 AMPANG, SELANGOR  
 TEL: 03-4296 5355 FAX: 03-4295 6005  
 E-MAIL: [sales@jatikhas.com](mailto:sales@jatikhas.com) WEBSITE: [www.jatikhas.com](http://www.jatikhas.com)

## ALUMINIUM TOP HOT PLATE



### Features

- ◆ Aluminium hot plate, uniform surface temperature
- ◆ Heating up fast and uniform, heating power is stepless adjustable

### Specifications

Model		BHP-6A	BHP-7A
	Working Time	Continuous	
	Max. Temp. (°C)	350	
Structure	Outer Shell	Iron plate, surface spraying	
	Hot Plate	Aluminum top	
	Heating Power (W)	1800	2000
Controller	Continuous Operation	Yes	
	Setting Method	Scale	
	Temperature Type	K	
	Exterior Size (W*L*H) (cm)	42*50*22	45*47*22
	Plate Size (mm)	350*450	400*400
	Plate Load	30 kg	30 kg
	Power Supply (50/60 Hz)	220 V	
	NW/GW (kg)	8.4/9.6	8.5/9.8