

BRANDON™

'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)

NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
68000 AMPANG, SELANGOR

TEL: 03-4296 5355 FAX: 03-4295 6005

E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

CERAMIC TOP HOT PLATE



Features

- ◆ Stretched shell and cover prevent leaking
- ◆ Heating and stirring can proceed simultaneously
- ◆ Enclosed hot plate with flame protection, fast warm-up and durability
- ◆ Heating power and stirring speed is stepless adjusted

Specifications

| Model | | BHP-4A |
|---------------|---|--------------------------------------|
| Function | Maximum Volume (mL) | 5000 |
| | Working Time | Continuous |
| | Maximum Temperature (°C) | 380 |
| Structure | Exterior | Iron plate welding, surface spraying |
| | Heating Plate | Ceramic Coated |
| | Heating Power (kW) | 0.6 |
| Controller | Continuous Operation | Yes |
| | Setting Mode | Scale |
| Specification | Exterior Size (W*L*H) (mm) | 320*190*125 |
| | Packing Size (W*L*H) (mm) | 390*240*230 |
| | Plate Size (mm) | 190*190 |
| | Plate Load | 10 kg |
| | Power Supply (50/60 Hz) Current rating | 110 V or 220 V |
| | NW/GW (kg) | 5/6 |