



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com



Safe storage of flammable and combustible liquids, in accordance with the requirements of the relevant regulations








Galvanized shelves for flammable safety cabinet. PP tray for corrosive safety cabinet.

FEATURES

- ◆ Safety cabinets are constructed of sturdy 18-gauge (1-mm) thick double-wall, welded steel with 1-1/2" (38-mm).
- ◆ Galvanized steel shelves direct spills to back and bottom of a leakproof 50mm/150mm bottom sump.
- ◆ Durable lead-free epoxy/polyester powder-coat finish provides excellent chemical resistance capacity.
- ◆ Three-point stainless steel bullet latching system makes it work safely every time.
- ◆ Equipped with Adjustable legs and Grounding wires

H x W x D (mm)	1650 x 1090 x 460
Capacity (Gal/Ltr)	45/170
Adjustable Shelves (pcs)	2
Shelves capacity (kg)	100
Door type	Double door
Model- Yellow	ZYC0045
Model- Blue	ZYC0045B

-  Yellow for flammable liquids
-  Red for paints, inks, and other combustible liquids
-  Blue for corrosive liquids
-  White for toxic substance
-  Green for pesticide