

BRANDON™

'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)

NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE, 68000 AMPANG, SELANGOR

TEL: 03-4296 5355 FAX: 03-4295 6005

E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

CONSTANT-TEMPERATURE SHAKING INCUBATOR



Features

- ◆ Large-screen LCD, shows temperature, speed, working time. Menu operation interface, easy to operate.
- ◆ Microcomputer control temperature and oscillation frequency with timing functions, Built-in power-off protection can automatically resume operation after power on.
- ◆ Unique air duct breeze circulation, Large observation window.
- ◆ Speed control : Intelligent Brushless DC motor, High speed accuracy, Slow start design, Auto lock if speed out of control. Adjustable starting speed to prevent the splash of shake flask liquid.

Specifications

MODEL		BSI-50B	BSI-70B
Shaking Mode		Cyclotron Oscillation	
Function	Amplitude	Φ20mm	
	Temp. Resolution	0.1 °C	
	Temp. Motion	±0.1°C	
	Temp. Uniformity	±1°C	
	Temp. Range	Ambient +5 ~ 65°C	4 ~ 65°C
	Rotary Speed	30 ~ 300rpm	30 ~ 300rpm
	Rotary Speed Precision	±1rpm	
Structure	Chamber Material	Mirror Stainless Steel	
	Outer Shell	ABS	
	Observation Window Material	Imported resin	
	Thermal Insulation Material	Polyurethane	
	Cooling System	No	Yes
	Coolant	No	R134a
	Heater	Stainless steel heating tube	
	Drive Mode	Three-axis level	
	Oscillating Plate	Spring Wire Rack	
	Air Circulation Mode	Horizontal convection	
	Power Rating	0.6kW	0.8kW



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

Controller	Temp. Control Mode	LCD PID intelligent control	
	Cyclotron Control Mode	Automatic changes of rotation direction	
	Setting Mode	Touch button setting	
	Temp. Display Mode	LCD display	
	Speed Display Mode		
	Timer	0 ~ 9999min	
	Operation Function	Fixed value operation, timing function, bidirectional automatic switchover, auto stop, quick stop start	
	Temperature Sensor	PT100	
	Speed Sensor	Hoare	
	Gating System	Yes	
Additional Function	Deviation correction, menu keys locked, power-off memory		
Safety Device		over temperature alarm, menu locked, gating system, power-off protection, plugging turn protection	
Specification	Inner Chamber Size (W*L*H)(mm)	455*380*290	560*390*320
	Exterior Size (W*L*H)(mm)	500*640*480	600*770*500
	Packing Size (W*L*H)(mm)	550*690*530	940*780*655
	Shelf Size	380*320	480*315
	Max capacity of Shelf for single specification (ml/pc)	100ml*16/250*10/500*6	100ml*28/250*15/500*8/1000*6
	Number of flask clamp as standard (ml/pc)	250*8	250*12
	Inside Volume	50L	70L
	Shelf Load Bearing	≤10kg/layer	≤15kg/layer
	Shelf Layer	1	
	Power Supply (50/60 Hz) Current Rating	AC220V/2.5A	AC220V/3.5A
	NW/GW (kg)	25/35	65/85