

BRANDON™

'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)

NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
68000 AMPANG, SELANGOR

TEL: 03-4296 5355 FAX: 03-4295 6005

E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

DESKTOP CONSTANT-TEMPERATURE DRYING OVEN/INCUBATOR



Features

- ◆ In the work room, the hot air is convective naturally and cyclically.
- ◆ Digital display microcomputer intelligent PID temperature controller.
- ◆ New anti-hot handle.

Specifications

MODEL		DRYING OVEN	INCUBATOR
		BHL-25AB	BP-25AB
Cycle Mode		Natural convection	
Function	Temp. Range	Ambient +10-300°C	Ambient +5-65°C
	Temp. Resolution Ratio	0.1 °C	
	Temp. Motion	±1 °C	±0.5°C
	Temp. Uniformity	±3.5%	±2.0°C
Structure	Inner Chamber	Stainless steel	
	Outer Shell	Cold rolling steel electrostatic spraying exterior	
	Insulation Layer	Compound silicate heat preservation board	
	Heater	Nickel chromium alloy heating wire	
	Power Rating	0.6kW	0.18kW
	Exhaust hole	φ28mm top (with function of test hole)	



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

Controller	Temp. Control Mode	Two temperature section intelligent PID	
	Temp. Setting Mode	Touch button setting	
	Temp. Display Mode	Measuring temperature: display on the four digital tube upper row; setting temperature: display on the lower row	
	Timer	0-9999 min (with timing wait function)	
	Operation Function	Fixed temperature operation, timing function, auto stop	
	Additional Function	Sensor deviation correction, temperature overshoot self-tuning, internal parameter locking, power-off parameter memory	
	Sensor	K model	CU50
Safety Device		Over temperature sound-light alarm	
Specification	Inner Chamber Size (W*L*H)(mm)	260*260*260	
	Exterior Size (W*L*H)(mm)	380*455*568	
	Packing Size (W*L*H)(mm)	500*530*657	
	Volume	18L	
	Shelf Number	6	
	Load per Rack	15kg	
	Shelf Space	29mm	
	Power Supply (50/60Hz)	AC220V/2.7A	AC220V/0.8A
	NW/GW (kg)	27/35	
Accessory	Shelf	2	
Optional Accessories		Shelf	