



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)  
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,  
 68000 AMPANG, SELANGOR  
 TEL: 03-4296 5355 FAX: 03-4295 6005  
 E-MAIL: [sales@jatikhas.com](mailto:sales@jatikhas.com) WEBSITE: [www.jatikhas.com](http://www.jatikhas.com)

## DIGITAL DISPLAY HEATING MANTLE



## Features

- ◆ PID Single Screen Temperature controller. Fast heating, temperature uniformity and high precise temperature control
- ◆ Several rows type, produce as customer inquiry, beautiful shape, easy to operate and effective
- ◆ C/w PT100 Temperature sensor, with adjustable sensor rack and stainless rods

## Specifications

Model		BHM-II-D				
Function	Capacity (mL)	250	500	1000	2000	5000
	Heating Mode	Surface heat-conduction				
	Surface Max. Temp. (°C)	450				
	Temperature Control Range	RT +20 ~ 250 °C				
Structure	Heating Element	Nichrome wire				
	Shell	Cold rolling sheet stretching and spraying exterior				
	Heat Power (kW)	0.15	0.25	0.35	0.45	0.80
Controller	Working Mode	Continuous				
Specification	Exterior Size (mm)	220*165	280*220	330*230	350*250	
	Packing Size (mm)	230*230*195	290*290*300	340*340*350	390*390*310	
	Voltage (V)	100 ~ 120 or 200 ~ 240				200 ~ 240
	N.W. (kg)	2.5	5.5	6.5	8	
	G.W. (kg)	2.8	6	7	9.2	