

ELECTRONIC CONTROL HEATING MANTLE



Features

- ◆ Electronic temperature regulation mode with characteristics of stepless heat adjustable knob.

Specifications

Model		BHM-I-A						
Function	Capacity (mL)	100	250	500	1000	2000	5000	
	Heating Mode	Surface heat-conduction						
	Surface Max. Temp. (°C)	450						
Structure	Heating Element	Nichrome wire						
	Shell	Cold rolling sheet stretching and spraying exterior						
	Heat Power (kW)	0.10	0.15	0.25	0.35	0.45	0.80	
Controller	Working Mode	Continuous						
Specification	Exterior Size (mm)	220*165		280*220	330*230	350*250		
	Packing Size (mm)	230*230*195		290*290*300	340*340*350	390*390*310		
	Voltage (V)	100 ~ 120 or 200 ~ 240					200 ~ 240	
	N.W. (kg)	2.5		5.5	6.5	8		
	G.W. (kg)	2.8		6	7	9.2		