

List of BRANDON Products

- Hotplate
- Hotplate Stirrer
- Overhead Stirrer
- Heating Mantle
- Water Baths
- Water Baths Shaker
- Orbital Shaker
- Shaking Incubator
- Incubator
- Drying Oven
- Muffle Furnace
- Chemical Safety Cabinet
- Laminar Flow
- Lab Benches
- Petri Dish
- Laboratory Glassware
- Stainless Steel Custom-Made/Custom Designed Products



Jatikhas Sdn Bhd

No. 10 & 10-1, Jalan Avenue 1, Avenue Ampang,
68000 Ampang, Selangor.

Tel: 03-4296 5355

Fax: 03-4295 6005

sales@jatikhas.com

jatikhas.sales@gmail.com

PIC: Vernon Jeremy Goonting
Anna Lauw

www.jatikhas.com

BRANDON™


General Laboratory Equipment and
Laboratory Design

**“High Quality Products
at Economical Price”**

~ 18 Months Warranty ~



JATIKHAS SDN BHD



LABORATORY SUPPLIES

BRANDON Equipment

HOTPLATE & HOTPLATE STIRRER



Aluminium Top Hotplate
BHP-6A/7A
• Plate Size:
350*450 (6A)
400*400 (7A)
• Max. Temp.: 350°C



Ceramic Top Hotplate
BHP-4A
• Max. Temp.: 380°C
• Plate Size:
190*190mm



Analog Ceramic Hotplate Stirrer
BHS-4
• Max. Temp.: 380°C
• Speed:
200-2000rpm
• Plate Size:
190*190mm



Digital Ceramic Hotplate Stirrer
BHS-II-4D
• Max. Temp.: 350°C
• Speed:
200-2000rpm
• Plate Size:
190*190mm

**comes with PT100 sensor

OVERHEAD STIRRER & HEATING MANTLE



Overhead Stirrer
BOHS-II-6D
• Capacity: 20L
• Speed:
100-2500rpm



Digital Display Heating Mantle
BHM-II-D
• Capacity: 250-5000mL
• Max. Temp.: 450°C



Electronic Control Heating Mantle
BHM-I-A
• Capacity: 100-5000mL
• Max. Temp.: 450°C



Extraction Mantle
BMSH-98-IV-B / BMSH-98-V-B
• With/Without Stirring Mode
• Capacity: 250-1000mL (4/6 holes)
• Max. Temp.: 450°C

**comes with PT100 sensor

WATER BATHS & SHAKER



Transparent Water Baths
BTWB-98-IV
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$
• Temp. Range:
RT+5-100°C



Thermostatic Water Baths
BWB.420A/600A (2/6 holes)
• With/Without Test Tube Rack & Cover Circle (AI/AII Model)
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$
• Temp. Range: RT+5-100°C



Thermostatic Water Baths (Stretching Inner Chamber)
BWB-2000-IIIIL (1-6 holes)
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$
• Temp. Range: RT+5-100°C



Constant Temperature Shaking Incubator
BSI-50B/70B
• With/Without Cooling System
• Rotary Speed Range: 30-3000rpm
• Temp. Accuracy: $\pm 0.1^\circ\text{C}$
• Temp. Range: RT+5-65°C
• Timer: 0-9999min



Multifunction Water Baths Shaker
BWS-30
• Rotary Speed Range: 20-2000rpm
• Temp. Accuracy: $\pm 0.1^\circ\text{C}$
• Temp. Range: RT+5-100°C
• Power Rating: 1.5kW
• Timer: 0-9999min



Orbital Shaker
BDS-O20
• Rotary Speed Range: 30-2000rpm
• Timer: 0-9999min (with timing function)

INCUBATOR



High Precision Constant-Temperature Incubator
BCI-II
• Capacity: 45-210L
• Temp. Range: RT+5-65°C
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$



Constant-Temperature Incubator
BCI-I-30BE/45BE/65BE/85BE/125BE/230BE
• Capacity: 30-230L
• Temp. Range: RT+5-120°C
• Temp. Accuracy: $\pm 1^\circ\text{C}$



Desktop Constant-Temperature Incubator
BP-25AB
• Capacity: 18L
• Temp. Range: RT+5-65°C
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$



Cooling Incubator
BCLI-70BIII/150BIII/250BIII
• Temp. Range: 0-65°C
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$
• Power: 300-800W
• Max. Load per Rack: 15kg



Cooling Incubator
BCLI-70BIV/150BIV/250BIV
• Temp. Range: 0-65°C
• Temp. Accuracy: $\pm 0.5^\circ\text{C}$ up to $\pm 1^\circ\text{C}$
• Power: 800-1200W
• Max. Load per Rack: 15kg
• Optional series with LCD display & sensor fault alarm (BL Model)

DRYING OVEN & MUFFLE FURNACE



Desktop Constant-Temperature Drying Oven
BHL-25AB
• Capacity: 18L
• Temp. Range: RT+10-300°C
• Temp. Accuracy: $\pm 1^\circ\text{C}$



Forced Air Drying Oven
BDO-II-30BE/45BE/65BE/125BE/230BE
• Capacity: 30-230L
• Temp. Range: RT+10-300°C
• Temp. Accuracy: $\pm 1^\circ\text{C}$



Big Forced Air Drying Oven
BDO-4AB/5AB/6AB/7AB/8AB
• Capacity: 640-3070L
• Temp. Range: RT+10-250°C
• Temp. Accuracy: $\pm 1^\circ\text{C}$



Ceramic Fiber Muffle Furnace
BMF3-2.5-12/5-12/8-12/10-12
• Capacity: 2-16L
• Temp. Range: 100-1200°C
• Temp. Accuracy: $\pm 1^\circ\text{C}$

CHEMICAL SAFETY CABINET & LAMINAR FLOW

CE EN 14470-1 OSHA EPA SRC COMPLIANT



Chemical Safety Cabinet
ZYC0045/0060/0030B/0045B
• Capacity: 30/114-60/227 gal/L
• 2 pcs of Adjustable Shelves
• Shelves Capacity: 100kg
• Equipped with adjustable legs & grounding wires



Laminar Flow Cabinet
BLC-1D/1S/2D/2S
• Airflow: Vertical
• Collection Efficiency: $\geq 99.9\%$ of $0.5\mu\text{m}$ particle
• Wind Speed Range: 0.2-0.6m/s

Galvanized shelves for flammable safety cabinet, PP tray for corrosive safety cabinet.

Lab Benches



- Laboratory benches, furniture fitting
- LEV testing and DOE submission for fume cupboard, alsident arm and canopy

Petri Dish



- Made of clear polystyrene, it allows for disposability and durability, while keeping clarity.
- Easy to stock & grip, and have anti-slip ring which offers stability.

Laboratory Glassware



Stainless Steel Custom-Made

