

## OVERHEAD STIRRER



### Features

- ◆ Brushless DC motor drive, noise and maintenance free
- ◆ Stainless steel head restraints, more durable. Stirrer with digital type has three kinds of rotation mode as positive rotation, inversion, reciprocal and mutual conversion

### Specifications

| Model         |  | BOHS-II-6D    |
|---------------|--|---------------|
| Function      | Max. Volume (L)                        | 20            |
|               | Digital                                | LCD           |
|               | Speed (r.p.m.)                         | 100 ~ 2500    |
|               | Working Time                           | Continuous    |
|               | Minimum Revolution (r.p.m.)            | 100           |
| Structure     | Exterior                               | Plastic shell |
|               | Driving Mode                           | Motor         |
|               | Motor Power (kW)                       | 0.06          |
| Controller    | Continuous Operation                   | Yes           |
|               | Setting Mode                           | Digital       |
|               | Speed Display Mode                     | Digital       |
|               | Timer                                  | Yes           |
|               | Operation Function                     | Yes           |
| Specification | Exterior Size (W*L*H) (mm)             | 115*95*240    |
|               | Packing Size (W*L*H) (mm)              | 400*250*210   |
|               | Plate Size (mm)                        | 200*300       |
|               | Plate Load                             | 50 kg         |
|               | Power Supply (50/60 Hz) Current Rating | 100 ~ 240 V   |
|               | NW/GW (kg)                             | 6/7           |