



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

OVERHEAD STIRRER



Features

- ◆ Brushless DC motor drive, noise and maintenance free
- ◆ Stainless steel head restraints, more durable. Stirrer with digital type has three kinds of rotation mode as positive rotation, inversion, reciprocal and mutual conversion

Specifications

Model		BOHS-II-6D
Function	Max. Volume (L)	20
	Digital	LCD
	Speed (r.p.m.)	100 ~ 2500
	Working Time	Continuous
	Minimum Revolution (r.p.m.)	100
Structure	Exterior	Plastic shell
	Driving Mode	Motor
	Motor Power (kW)	0.06
Controller	Continuous Operation	Yes
	Setting Mode	Digital
	Speed Display Mode	Digital
	Timer	Yes
	Operation Function	Yes
Specification	Exterior Size (W*L*H) (mm)	115*95*240
	Packing Size (W*L*H) (mm)	400*250*210
	Plate Size (mm)	200*300
	Plate Load	50 kg
	Power Supply (50/60 Hz) Current Rating	100 ~ 240 V
	NW/GW (kg)	6/7