



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)  
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,  
 68000 AMPANG, SELANGOR  
 TEL: 03-4296 5355 FAX: 03-4295 6005  
 E-MAIL: [sales@jatikhas.com](mailto:sales@jatikhas.com) WEBSITE: [www.jatikhas.com](http://www.jatikhas.com)

## SEVERAL ROWS HEATING MANTLE (WITH/WITHOUT STIRRING)



## Features

- ◆ BMSH-98-IV-B: Electronic temperature regulation heating mode.
- ◆ BMSH-98-V-B: Electronic temperature regulation heating mode. Heating and stirring can proceed simultaneously; stepless speed adjustment.

## Specifications

Model		BMSH-98-IV-B/BMSH-98-V-B											
Function	Capacity (mL)	2X 100	2X 250	2X 500	2X 1000	4X 100	4X 250	4X 500	4X 1000	6X 100	6X 250	6X500	6X1000
	Heating Mode	Surface heat-conduction											
	Surface Max. Temp.	450°C											
	Temp. Range	Ambient+vf20 ~ 250°C											
	Speed Range	BMSH-98-IV-B: No						BMSH-98-V-B: 0~1600rpm					
Structure	Heating Element	Nickel Chromium alloy wire											
	Stir Bar	BMSH-98-IV-B: No						BMSH-98-V-B: PTFE					
	Outer Shell	Cold rolled steel plate, surface spraying											
	Motor Power	BMSH-98-IV-B: No						BMSH-98-V-B: 0.01kW					
	Heating Power (kW)	2X0.1	2X 0.15	2X 0.25	2X 0.35	4X 0.1	4X 0.15	4X 0.25	4X 0.35	6X 0.1	6X 0.15	6X 0.25	6X 0.35
Controller	Working Mode	Continuous											
Specification	Exterior Size (mm)	280* 140* 106	320* 160* 170	360* 180* 180	420* 210* 200	530* 140* 160	610* 160* 170	690* 180* 180	810* 210* 200	820* 140* 160	940* 160* 170	1060* 180* 180	1240* 210* 200
	Packing Size (mm)	320* 180* 190	360* 190* 200	400* 220* 210	570* 180* 190	570* 180* 190	650* 190* 200	730* 220* 210	850* 250* 230	860* 180* 190	980* 190* 200	1100* 220* 210	1280* 250* 230
	Voltage	100 ~ 120V or 200 ~ 240V						200 ~ 240V					
	N.W. (kg)	4	4	4	7	8	8	8	14	12	12	12	21
	G.W. (kg)	5	5	5	8	10	10	10	16	15	15	15	24