



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

THERMOSTATIC WATER BATHS



Features

- ◆ **BWB.420A** model cover for the two holes water bath, can be used as a digital two holes water bath.
- ◆ **BWB.600A** model cover for the six holes water bath, can be used as a digital six holes water bath.

Specifications

Model		BWB.420AI	BWB.420AII	BWB.600AI	BWB.600AII
		1 Row & 2 Holes		2 Rows & 6 Holes	
Function	Temp. Range	Ambient+5~65°C	Ambient+5~100°C	Ambient+5~65°C	Ambient+5~100°C
	Temp. Resolution Ratio	0.1°C			
	Temp. Accuracy	±0.5°C			
	Temp. Uniformity	±1.0°C			
Structure	Inner Chamber	High quality hairline stainless steel plate			
	Outer Shell	Cold rolling steel electrostatic spraying exterior			
	Cover	High quality hairline stainless steel plate			
	Heater	Stainless steel heating tube			
	Cover Circle	Without	With	Without	With
	Test Tube Rack				
	Power Rating (kW)	0.5	1.0	0.75	2.0
Controller	Temp. Control Mode	PID			
	Temp. Setting Mode	Touch button setting			
	Temp. Display Mode	Measuring temperature: display on the three digital tubes upper row; Setting temperature: display on the lower row			
	Timer	0~9999min (with timing wait function)			
	Operation Function	Fixed temperature operation, timing function, auto stop.			
	Program Function	No			
	Additional Function	Deviation Correction, menu key lock, power-off compensation, power failure memory			
	Sensor	CU50			
Safety Device		Over temperature alarm			
Specification	Inner Chamber size (W*L*H) (mm)	420*180*140		600*300*200	
	Exterior Size (W*L*H) (mm)	570*230*270		750*350*360	
	Packing Size (W*L*H) (mm)	670*320*335		860*460*380	
	Volume (L)	11		36	
	Power Supply (50/60Hz) Current rating	AC220V/2.3A	AC220V/4.5A	AC220V/6.8A	AC220V/9.1A
	NW/GW (kg)	8/9.5	8.5/10	17/19	18/20