



'One Stop Centre' supplier for all scientific equipment, glassware and safety equipment for local and international industries

JATIKHAS SDN BHD (402338-A)
 NO. 10 & 10-1, JALAN AVENUE 1 AMPANG AVENUE,
 68000 AMPANG, SELANGOR
 TEL: 03-4296 5355 FAX: 03-4295 6005
 E-MAIL: sales@jatikhas.com WEBSITE: www.jatikhas.com

TRANSPARENT WATER BATHS



Features

- ◆ Using innovative design structure, effectively prevent steam leakage.
- ◆ Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
- ◆ Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.

Specifications

MODEL	BTWB-98-IV
Voltage	200-240V/100-120V
Power (kW)	0.8
Temp. Accuracy (°C)	±0.5
Temp. Range (°C)	Ambient+5~100
Temp. Sensitivity (°C)	≤±1
Shelves Height	Distance from bottom 100mm
Shelves Holes Diameter (mm)	∅ 5*104
Inner Chamber Size (L*W*H) (mm)	280*220*150
Exterior Size (L*W*H) (mm)	400*260*190
Packing Size (L*W*H) (mm)	450*350*270
N.W. (kg)	3.5
G.W. (kg)	5